HYDROCHLOR WHIRLPOOL ANTISEPTIC



Manufacturer of a Comprehensive Line of PROFESSIONAL HEALTHCARE & THERAPY EQUIPMENT

AVAILABLE IN EASY TO USE WATER DISSOLVABLE PACKETS

Our HydroChlor water dissolvable packets are highly effective with hydrotherapy treatments. Simply fill the whirlpool to desired level, drop in the required size of HydroChlor[™] packet, agitate 2-3 minutes. It's that easy!

HYDROCHLOR

- Contains 100% Chloromine-T, proven to be an effective way to sanitize.
- Is available in pre-measured water-dissolvable packets or in bottles. Please see below for specific sizes and model numbers.

Model	#	Description

HC-5	5 gram, 8 units per box
HC-15	15 gram, 48 units per box
HC-30	30 gram, 36 units per box
HC-50	50 gram, 24 units per box
HC-80	80 gram, 24 units per box
HC-200	200 gram HydroChlor bottle
HC-320	320 gram HydroChlor bottle (sold in case qty. 12 to a case)
*HC-5LB	5 lb HydroChlor bottle (sold in case qty. 4 to a case)

- *HC-5LB-1 5 lb HydroChlor bottle (sold individually)
- * Non-returnable



How much HydroChlor do I need?

- If your tank holds 3-5 gallons \longrightarrow Use a 5 gram packet
- If your tank holds 26-40 gallons —
- If your tank holds 44-66 gallons \longrightarrow Use a 50 gram packet

- If your tank holds 282-423 gallons \longrightarrow Use a 320 gram bottle

- \longrightarrow Use a 30 gram packet
- If your tank holds 176-264 gallons Use a 200 gram bottle





Dry Heat





190. Boul. St-Elzéar Ouest Laval (Québec) H7L 3N3 T 450 669.8985 | 1 800 361.0877 F 450 669.9532 1888 935.7001 www.medicaltronik.ca

Paraffin Cold Therapy Splint Pans

P.O. Box 3527 • City of Industry, CA 91744-0527, U.S.A. • (800) 782-7706 • (626) 968-6681 • www.whitehallmfg.com



