



# ELECTRIC TREATMENT TABLE

# CONFYCARE

## MASSAGE

- Electric height adjustable by footbar system (19"- 36.8")
- Tilting headpiece : -75° to 30°
- 4 multidirectional casters
- 7" extendable/removable extension
- Cushions thickness: 2.8"
- Loading weight: 495 lb (225 kg)
- Certification: tested according to the Canadian standard by QPS  
CE, FDA, ISO 13485-2016, EN60601
- Warranty:
  - Motor/Control box: 5 years
  - Frame: 5 years
  - Cushions: 3 years
  - Parts/Cylinder: 1 year
  - Labor: 1 year



### AVAILABLE COLOURS



BLACK



CHARCOAL



NAVY  
BLUE

DISTRIBUTED BY



190, Boul. St-Elzéar Ouest  
Laval (Québec) H7L 3N3  
T 450 669.8985 | 1 800 361.0877  
F 450 669.9532 | 1 888 935.7001

[www.mtm.ca](http://www.mtm.ca)