



ELECTRIC CHIROPRACTIC TABLE



- Electric elevation (17" - 31")
- Adjustable width tilting headpieces: 1.25"
- Tilting headpiece: -20° to +25 °
- Raised headpiece: 4.75"
- Cervical drop
- Thoracic drop
- Lumbar drop
- Pelvic drop
- Raised pelvic
- Cushions thickness: 2,4"
- Can be moved with retractable wheel system
- Loading weight: 495 lb (225 kg)
- Certification:
Tested according to the Canadian standard by QPS, CE, FDA, ISO9001-2008, ISO13485-2007
- Warranty:
 - Motor/Control box: 4 years
 - Frame: 4 years
 - Cushions: 2 years



AVAILABLE COLOURS

- | | | | |
|-----------------|-------|----------|-----------|
| CUSHIONS | BLACK | CHARCOAL | NAVY BLUE |
| FRAME | IVORY | GRAPHITE | |

DISTRIBUTED BY



190, Boul. St-Elzéar Ouest
Laval (Québec) H7L 3N3
T 450 669.8985 | 1 800 361.0877
F 450 669.9532 | 1 888 935.7001

www.mtm.ca