

190, Boul. St-Elzéar Ouest, Laval (Québec) H7L 3N3 T 450 669.8985 | 1 800 361.0877 www.mtm.ca

## Cranberry®



CR3230 Series 200 gloves/box 10 boxes/case XS: CR3235
S: CR3236
M: CR3237
L: CR3238
XL: CR3239

Gloves

CARBON®

Nitrile Powder Free Examination Gloves

- Enhanced fingertip texture for superior handling and increased control
- Distinct black color to minimize staining during use
- 200 space saver pack reducing storage space and packaging waste

Product Information					
Material	Nitrile		Finger Thickness	0.10 mm ± 0.02 / 3.4 mil	
Color	Black		Palm Thickness	0.07 mm ± 0.02 / 3.2 mil	
Length (mm)	240				
Texture	Fingertip		Powder Content	Powder-Free	
Sterile	No		Tensile Strength	14 MPa	
Product Sizing					
Palm Width (mm)	XS (CR3235)	Small (CR3236)	Medium (CR3237)	Large (CR3238)	XL (CR3239)
	73 - 79	81 - 87	91-97	102 -108	110 - 116

## Ecological Information

Bio-Compatibility: The chemical formulation of the gloves and surface lubricating materials do not contain any substances normally known to be harmful to the user or to any person with whom the gloves contact.

Other information

To prolong product life, gloves should be stored in a cool, dry location, located away from heat and light.

The data in this Material Safety Data Sheet relates to the specific material designated herein and does not relate to the use in combination with other material or in any process.