NeuroTrac[®] MyoPlus 2 and 4 Pro

The Professional Choice

MyoPlus 2 PRO

Dual channel STIM / EMG / ETS

Ideal for professional use in clinic for pelvic floor disorders. Assessment and training tool.




MyoPlus 4 PRO

Four channel STIM and Dual channel EMG / ETS Ideal for professional use in clinic for muscular rehabilitation. Assessment and training tool.



- * The Professional EMG tool for Biofeedback training and neuro-muscular assessment.
- Advanced ETS used for neuro-reeducation: the graph shows the muscular effort and the stimulation is triggered automatically to support the effort.
- Pre-set and customised programmes for Pelvic Floor and Incontinence, Muscular Rehabilitation and Sports.
- * Upto 10 metres of wireless (Bluetooth) connection for NeuroTrac Software





DISTRIBUTED BY



