



THE GO EVERYWHERE, TREAT EVERYTHING LASER THERAPY SOLUTION.

 $\bigcirc \bigcirc$

07:00

A 251

START TREATMENT

×251 (+ 251)

200 J

(25%

Total Joules

1000 J

 \mathcal{O}

EDIT

With a lightweight yet sturdy design, the H1 delivers powerful results and simple operation. It's built on our advanced Android platform so that you can easily and quickly select an expert treatment protocol that will get your patients healing better, faster.

> 05:30 PHASE 3 02:32 8 W 20.000 Hz Body Part Treating



15 Watts in CW (16 Watt ISP Peak) More power than most devices twice its size.



165 Expert Treatment Protocols

Designed and tested by Veterinarians for Veterinarians. Plus, create and save your own custom protocols.





Portable & Sturdy

Lightweight yet durable construction allows you to take the laser everywhere you need it.



Ease of Use


Intuitive touchscreen interface allows you to quickly set-up an expert treatment protocol.



Best-In-Class Warranty Lifetime diode warranty and 3-year manufacturer's warranty.

3 WAVELENGTHS

Simultaneous or independent application effectively penetrates tissue for excellent results.





SPECIFICATIONS

Wavelengths	3 Wavelengths 650nm I 810nm I 915nm
Power	15 watts in CW
Peak Power	16 watts in ISP
Hand Piece Lens	25 mm zoom - Open lens
Pulsing	1 to 20,000 Hz & ISP
Beam Profile	Homogeneous beam
Diode	High powered solid-state diode module
Battery	 • Up to 1 hours (in-use) • 30-minute recharge time • 7 hours (stationary) • Interchangeable battery
Technology Platform	 7 Inch 1920x1200 HD touch screen Android platform Wi-Fi enabled software updates
Weight	4.4 lbs.
Dimensions	8.9"(L) X 5.9"(W) X 5.9"(H)
Warranty	Lifetime diode warranty and 3-year manufacturer's warranty



