

Mobile 2 RPW



Mobile 2 RPW



Chattanooga's Radial Pressure Wave technology is mobile.

The Mobile 2 RPW is compact and light enough to be easily carried while also delivering a host of features and benefits.

When there is need for a powerful yet portable Radial Pressure Wave device in the clinic or on the sports field, the Chattanooga Mobile 2 RPW is the ideal choice. A lightweight and compact design together with a quiet built-in compressor make it perfect for portable use. A user-friendly experience is at the heart of the Mobile 2 RPW's design. The simple, ergonomically laid out controls enable all of the key parameters to be easily changed, while a clear, illuminated LED display removes ambiguity when monitoring pressure, frequency and number of pulses. Pulses can be initiated via a trigger mechanism on the handpiece.



User Profile

This Radial Pressure Wave clinical unit is ideally suited to the following applications:

- Orthopaedics
- Rheumatology
- Physical Medicine
- Physiotherapy
- Sport Clubs
- Acupuncture
- Aesthetic
- Veterinary
- Podiatrists
- Chiropractors

Indications

- Myofascial pain
- Muscle tendon tissue activation
- Tendon attachment & site illness
- Acupuncture Shockwave Therapy

The Mobile RPW D-Actor Applicator allows quick and efficient uses of applicator transmitters for desired therapeutic effects during treatment.

D-ACTOR® Applicator

Applications:

- Myofascial trigger points
- Other disorders of tendon insertions
- Activation of muscle and connective tissue
- Acupuncture shock wave therapy
- Cellulite treatment (stages I – III)
- Muscle and connective tissue tightening
- Elimination of lymphatic congestion
- Body shaping
- Regulation of muscle tone
- Medical massage



Facts and figures

- Weight of the unit: Approx. 9 kg
- Dimensions (L x W x H): 32 x 32 x 21 cm
- Built in compressor
- Single and Continuous shocks
 - Frequency D-actor: 1 - 21 Hz
 - Frequency V-actor: 1 - 31 Hz
- Application pressure: 1 - 5 Bar
 - Highly customisable with precise adjustments of 0.1 increments.
- Oscillating D-Actor
- V-Actor Hand piece option
- USB connectivity



User Interface Features

1. Pressure screen 1 - 5 Bar
2. Pulse counter screen
3. Frequency screen
4. Intensity adjustment
5. Pulse counter reset
6. Frequency adjustment

Hand piece docking port



Ultrasound gel holder



D-ACTOR® Transmitters

C 15	Ro40	D 20	D 35	DI 15	F 15
0.45 mj/mm ²	0.38 mj/mm ²	0.48 mj/mm ²	0.46 mj/mm ²	0.63 mj/mm ²	0.16 mj/mm ²
0 - 40 mm	0 - 40 mm	0 - 50 mm	0 - 50 mm	0 - 60 mm	0 - 20 mm



Mobile 2 RPW

Ordering Information

Part Number	Description
2910	Mobile 2 RPW

Standard Accessories (included with Mobile 2 RPW as default)

2905	Chattanooga Mobile 2 RPW Unit Only
32000	R-SW Hand Piece Applicator Set
28178	Ro 40 'Beam' transmitter
29724	D20-S D-Actor * 20mm Transmitter
13-5182	Conductor™ Transmission Gel 250 ml
13-00060	CD User Manual

Optional Accessories

28701	Revision Kit
28178	Ro 40 'Beam' transmitter
31800	9 Transmitter Set In Case, Accessory
29801	A6 Point, Transmitter
29802	T10 Finger, Transmitter
29729	C15, 15 mm Cerama-X Transmitter, Black
29728	D115, 15 mm, Deep Impact Transmitter, Black
29726	F15 Focus-Lens 15 mm Transmitter
29724	D20-S D-Actor * 20mm Transmitter
28736	D20-T Transmitter
29725	D35-S, 35 mm, D-Actor Transmitter, Black
22375	Handpiece Cushion, Mobile 2 RPW
22376	Shaft Cushion, Mobile 2 RPW
28730	V-Actor Set
28740	V-Actor V25 Transmitter
28741	V-Actor V35 Transmitter
28745	Carrying Bag



MOTION IS MEDICINE®

DISTRIBUTED BY



190, Boul. St-Elzéar Ouest
Laval (Québec) H7L 3N3
T 450 669.8985 | 1 800 361.0877
F 450 669.9532 | 1 888 935.7001

www.mtm.ca