

## Computer Specifications for Commander Echo Downloader

The following computer system requirements must be met for JTECH's Commander Echo Downloader to function correctly. Although the software can run on systems that do not meet all the requirements, JTECH does not guarantee its performance.

	Minimum	Recommended
Operating System	32-bit or 64-bit Windows XP Professional, Windows Vista or Windows 7	32-bit or 64-bit Windows XP Professional, Windows Vista or Windows 7
CPU	1.0 Ghz	1.4Ghz+
Memory (RAM)	512MB	2GB+
Display	1024 x 768 px	1024 x 768 px
Hard Disk Space	500MB Free	1GB+ Free
Open Ports	1 x USB 2.0	1 x USB 2.0
Other	Microsoft Word® 2003 or newer	Microsoft Word® 2003 or newer

## **Important Notice:**

For customer service calls, it is necessary to have a telephone located within reach of the computer running the Commander Echo System. This is vital during installation and technical support calls as you will need to access the computer while working with our customer service staff. Without a telephone located near the computer, JTECH representatives and support staff will be unable to provide you with proper, timely service. If you have any further questions about computer system requirements for running Commander Echo Downloader, please contact your Sales Representative or JTECH Customer Service at (800) 985-8324.

Note: Windows 95, Windows 98, Windows 2000, Windows ME, Windows XP Home edition are not recommended or supported for Commander Echo Downloader.



