

Wireless Commander Echo Grip is the fastest, most versatile portable grip testing system for hand strength testing or voluntary maximal effort evaluations.



Advanced Features

- ▶ Performs 5-position, rapid exchange and sustained grip tests
- ▶ Displays maximum force, CV, right and left average force and bilateral deficit
- ▶ Evaluate voluntary maximal effort
- ▶ Establish baselines and monitor progress
- ▶ Adjustable test time with fatigue calculation for sustained grip tests
- ▶ Digital load cell technology provides greater reliability and durability than other common gauges
- ▶ Wireless console is independent of the testing device allowing the tester to review data during active testing without influencing the patient's level of effort



*Console required to operate device. Console and device sold separately.

DISTRIBUTED BY



190, Boul. St-Elzéar Ouest
Laval (Québec) H7L 3N3
T 450 669.8985 | 1 800 361.0877
F 450 669.9532 | 1 888 935.7001
www.medicaltronik.ca



JTECH MEDICAL
MEASURING HUMAN PERFORMANCE

© 2012 JTECH MEDICAL MM316_B

Technical Specifications

- Number of tests stored on Console:20
 - Number of repetitions per test: 1-8
 - Audible cues: yes
 - Load cell capacity: 200 lb.
 - Auto coefficient of variation calculations (CV): yes
 - Automatic bilateral deficit calculation: yes
 - Threshold settings: 3, 5, 10, 20 lb.
 - Units of measure (pounds default): . . .lb, kg, N
 - Automatic fatigue calculation: yes
 - Timed test for fatigue evaluation yes
 - Allows for specific test re-testing: Last rep only while still testing
 - Print options: Yes*
- * With purchase of Downloader Echo software



Product Highlights

- Console stores 20 bilateral tests with up to 8 repetitions per side
- Compatible with normative data developed by Mathiowetz and published in the American Medical Association (AMA) *Guides*
- Helps determine consistency of effort in Workers' Compensation and FCE exams
- Measures up to 200 pounds for very strong performers

COMMANDER
ECHO Downloader™

Optional Commander Echo Downloader (CM311) software enables you to create quick but meaningful reports from the data

you collected on your Commander Echo Console with the click of a button. The easy-to-use Commander Echo Downloader software interfaces directly with your Commander Echo Console eliminating dictation and transcription errors.



Other
Commander
Echo Devices*



*Some devices shown with accessories attached.

DISTRIBUTED BY



190, Boul. St-Elzéar Ouest
Laval (Québec) H7L 3N3
T 450 669.8985 | 1 800 361.0877
F 450 669.9532 | 1 888 935.7001
www.medicaltronik.ca



JTECI-MEDICAL
MEASURING HUMAN PERFORMANCE
© 2012 JTECH MEDICAL MM316_B