Wireless Muscle Tester



Quickly and accurately quantify muscle strength and identify areas of weakness with the wireless Commander Echo Muscle Tester.





Advanced Features

- Test protocols can be customized according to number of tests, starting force, Newtons, Kilograms or Pounds and number of repetitions
- ▶ Innovative pop-in "quick connect" system makes changing test pads fast and easy
- Quickly and accurately test muscle strength and identify areas of weakness with the tester and Commander Echo Console
- Quantify bilateral strength differences
- Identify patterns of weakness
- Provide objective measure of strength
- Wireless console is independent of the testing device allowing the tester to review data during active testing without influencing the patient's level of effort







*Console required to operate device. Console and device sold separately. Device shown with accessory attached.





Technical Specifications

* With purchase of Downloader Echo software

Included Accessories:

Flat Pad (8AC008) – The flat pad can be used for exerting force against flat items.

Curved Pad (8AC009) – The curved pad can be used for exerting force against items that are rounded.

Muscle Tester Cradle (9AC001) – The Muscle Tester Cradle provides a stable platform to rest the muscle tester.





Product Highlights

- Console records and stores up to 20 tests with up to 4 repetitions per side
- Stability during tests improved with "Axis Compensation" circuitry which guarantees 99% accuracy for off-center placement
- Reduces "roll-over" and tilting with a 2.11" high ergonomic design that fits in the palm
- Measures up to 150 pounds for very strong performers
- Automatically records peak force and calculates
 CV and bilateral strength deficits



Optional Commander Echo Downloader (CM311) software enables you

to create quick but meaningful reports from the data you collected on your Commander Echo Console with the click of a button. The easy-to-use Commander Echo Downloader software interfaces directly with your Commander Echo Console eliminating dictation and transcription errors.



Other











*Some devices shown with accessories attached.

DISTRIBUTED BY



190, Boul. St-Elzéar Ouest Laval (Ouébec) H7L 3N3 T 450 669.8985 | 1 800 361.0877 F 450 669.9532 | 1 888 935.7001 www.medicaltronik.ca

