





**VET TABLE**Elevating DR Table



**VET DENTAL** Intra-Oral Imaging System



**BOVINE & ZOO**DR Portable System



**EQUINE**DR Portable System



**SOFTWARE** iQ-View & Cloud System



# **EXOTIC & AVIAN | Hand-Held X-ray & Detector**

A lightweight and powerful battery operated, portable X-ray system that produces superior image quality. Designed for simplicity and ease of use, it is ideal for imaging exotic and avian animals.

### **Detector features**

- Removable handle
- Support Cloud system
- Wireless
- Auto image rotation
- Sleep mode
- Light weight durable design
- · Built-in image storage
- Shock sensor





# **MYVET TABLE | for Small Animals**

The MyVet table system is a revolutionary Veterinary Elevating X-Ray Table which delivers silent elevating to fulfill comfort and stability.

The MyVet table system is a cost effective, workflow efficient approach to digital radiographic imaging designed to improve patient workflow in veterinary practice of small animals.

- Completely Silent Elevating Table system
- Lowers to 10.5 inches off the floor
- Vertical 40 inches and 72 inches SID
- Horizontal Beam at 72 inches SID
- 17x17 Csl Flat Panel Detector
- Gives you the ability to do weight bearing exams
- Easy table operation by a foot switch

Charging Bucky

# **VET DENTAL | Intra-Oral Imaging System**

Our Intra-oral imaging system is the ideal veterinary digital imaging solution. It provides higher diagnostic value with superior image quality for precise and advanced intra-oral radiography.







Standard Equipment

# **BOVINE & ZOO | DR Portable System**

The Rayence portable package provides the ultimate in adaptability to its environment. The 1417 Series flat panel detector is the latest in innovative flat panel technology. The durable, ultra-thin and light weight design makes it the most ideal 14"x17" flat panel detector.

## **Standard Equipment**

- 1417 DR flat panel detector
- Laptop computer with pre-installed VetView software
- Portable carrying case with wheels
- Rechargeable 2-battery pack with a charger

## **Optional Equipment**

- Portable generator
- Tablet PC
- Mobile battery charger

# **EQUINE PACKAGE | DR Portable System**

The DR portable package for equine provides the ultimate in adaptability to its environment. It uses auto-trigger signal sensing technology that does not require any integration into an X-ray generator.

## **Standard Equipment**

- 1012 DR flat panel detector
- Laptop computer with pre-installed VetView software
- Portable carrying case with wheels
- Rechargeable 2-battery pack with a charger

### **Optional Equipment**

- Portable generator
- 14"x17" DR panel
- Tablet PC
- Mobile battery charger



(Optional)



Standard Equipment

# VET-SPECIFIC PLATFORM ZERO FOOTPRINT X-RAY WET-SPECIFIC PLATFORM X-RAY ULTRASOUND VETERINARY SERVICES VETERINARY PRACTICES, HOSPITALS AND GROUPS IMAGE ACCESS FROM ANYWHERE

# **SOFTWARE | iQ-VIEW & CLOUD**

iQ-View MyVet Edition is an easy-to use reading station, specialized for veterinarians. iQ-VIEW can easily be integrated with any practice management system (PMS), EMR or HIS. Even specialty tools for 3D post-processing, orthopedic templating, or stitching are closely interfaced with the viewer.





MyVet Cloud



## **Dealer Information**



