

Ordering Information

Part Number	Description
-------------	-------------

Primera TENS

77621	Primera TENS/NMES International
-------	---------------------------------

Complete kit includes manual, lead wires (2), garment belt, battery, and electrodes in the plastic case

Standard Accessories

77616	Primera Unit Sliding Back Cover
-------	---------------------------------

77622	User's Manual
-------	---------------

77623	Carrying Case
-------	---------------

77619	Leadwire Set - 2 each
-------	-----------------------

77620	Electrode Positioner Belt for the Back
-------	--

200001-001	Energizer 9-volt
------------	------------------

42192	Dura-Stick+ 5 cm Round Reorder Electrodes (Includes 10 packs of 4 electrodes in each pack)
-------	--

Technical Specifications

Number of Channels:	2
---------------------	---

Number of Programs:	14
---------------------	----

Constant Voltage:	Up to a resistance of 1000 ohms
-------------------	---------------------------------

Maximum Amplitude:	40mA @1000ohms, 250 µs
--------------------	------------------------

Output Intensity:	0-60, adjustable in 0.5 steps
-------------------	-------------------------------

Max Pulse Duration/width:	Up to 250 µseconds
---------------------------	--------------------

Max Frequency:	1,000 Hz
----------------	----------

Power Source:	PP3 Alkaline, 9V battery
---------------	--------------------------

Dimensions*:	108 x 65 x 23 mm
--------------	------------------

Weight:	70 g without battery, 100 g with battery
---------	--

Electrical Safety Class:	Class IIa, Type BF
--------------------------	--------------------

Safety Test:	IEC 60601-1, IEC 60601-1-2, IEC 60601-2-10
--------------	--

* Dimensions are listed as L x W x H

Waveforms



- Asymmetrical, rectangular bi-phasic with zero DC Current

Primera TENS™



Features

- Extensive range of programs – 8 TENS and 6 NMES Programs
- The HAN TENS Waveform*, further differentiates PRIMERA
- The PRIMERA features a proven circuit design with over 10 years of reliable performance
- Complete package includes: PRIMERA Unit, Carrying Case, Battery, Leadwires, Dura-Stick Electrodes, Back Support Belt and User's Manual
- 2 Year Warranty (Unit Only)

* Dr. Ji Sheng Han at the Neuroscience Research Institute at Beijing Medical University in Beijing China has done extensive research into electro-acupuncture and the role of endogenous opioid peptides in pain control for decades.

Dr. Han has demonstrated that his "HAN TENS Waveform" produce a synergistic analgesic effect for faster, greater and more long-lasting pain relief results.

