


maxill inc., St. Thomas, Ontario N5R 6C8 www.maxill.com 1-800-268-8633

# maxill VINYI

## **MEETS:**





Length 240 mm

### **ATTRIBUTES:**

- Powder free
- Premium quality
- Ambidextrous
- Short cuff: beaded
- Non-sterile
- Single use
- Latex free
- Smooth

# **DIMENSIONS:**

Box dimensions: 23.5 cm x 13.3 cm x 6.5 cm

Box count: 100 gloves per box Case pack: 10 boxes per case

SIZE	X-SMALL	SMALL	MEDIUM		LARGE	X-LARGE
Item Number	#51130	#51131	#51132		#51133	#51134
Width of Palm (mm)	75 ± 5	85 ± 5	95 ± 5		105 ± 5	115 ± 5
Total Length (mm)	240 min	240 min	240 min		240 min	240 min
Thickness: Palm (mm)	0.1 min	0.1 min	0.1 min		0.1 min	0.1 min
Thickness: Finger (mm)	0.1 min	0.1 min	0.1 min		0.1 min	0.1 min
PROPERTIES	BEFORE AGING			AFTER AGING		
Tensile Strength (MPa)	10.5 min			10.5 min		
Elongation (%)	350 min			350 min		
Force at Break (N)	4 min			4 min		



#### **Product of China**

The information provided is given in good faith and believed to be accurate, but no warranty or guarantee is given or implied. Conditions of use under which the product is being applied will vary from the conditions of normal laboratory efficacy testing. Therefore, each health care facility should conduct its own efficacy testing in order to verify the efficacy of the product under its own normal operating conditions to meet the specific requirements of the workplace.