

theratherm™

From the Makers of the *Original* Hydrocollator® Hot & Cold Products



Digital Moist Heating Pad

The Best Choice in Home Moist Heat Therapy

From the creators of the Hydrocollator® Moist Heat HotPac™, Chattanooga Group is pleased to introduce Theratherm, a digital electric moist heating pad designed for the convenience of home use. The digital moist heating pad provides temporary relief from pain associated with muscle spasm and inflammation caused by strain and tension.

The Theratherm Digital Moist Heating Pad is easier to use than the other electric heating pads on the market due to its unique digital hand control. Most competitive units require continuous contact with a spring-loaded switch to generate and maintain heat from the pad. The digital controller allows one to set the temperature and time as well as features a lockout mode to guard against remiss touching of the controls during treatment.

- **Digital hand control** allows operator to program, as well as monitor, the temperature and treatment time.
- **Lockout mode** protects against inconsistent changing of the control settings throughout treatment.
- **Allows patient to relax** and enjoy therapy without the hassle of maintaining contact with a hand held switch.
- **Therapeutic treatment temperature** between 88° and 166°F.
- **Treatment time** between 1 to 60 minutes.
- **Specially designed 100% cotton flannel cover.**

Available in three convenient sizes:

1032 – Standard

14" x 27" (36 cm x 68 cm)

1031 – Medium

14" x 14" (36 cm x 36 cm)

1030 – Small

7" x 15" (18 cm x 38 cm)

1035 – Retail

Floor Display

Includes Six Theratherms (two of each size)

